

Integrated Bioparticle Concentration for Improved Detection in Autonomous Systems

Sample Prep 2011, San Diego, CA



Andrew Page¹, David Alburty¹, Pam Murowchick², Zachary Packingham¹

¹InnovaPrep, LLC

²AlburtyLab, Inc.


Outline

- I. InnovaPrep/Hydrosol Concentration Background
- II. Performance
- III. Integration with the Evogen Dryclone
- IV. Integration for DHS AESaP
- V. Integration with MFSI M-BAND

InnovaPrep/Hydrosol Concentration Background




Automated Detection of Aerosolized Biological Weapons – **Standard Practice**



***75% to more than 99% of the sample is not analyzed!**

Automated Detection of Aerosolized Biological Weapons – **With Concentration**



***Delivery of nearly 100% of the collected signature is possible – providing orders of magnitude improvement in DL**


Implementation of a Hydrosol Concentrator

Advantages

- Orders of magnitude improvement in detection limit
- Increased freedom in selection of collector
- Increased freedom in selection of detector

Disadvantages


- It adds another component to the system
- It takes additional time
- Co-concentration of Inhibitors is possible



The InnovaPrep Concentration Process

- “Dead-end” filtration through a hollow fiber or flat membrane filter
- Particles larger than the filter pore size are retained while the permeate fluid passes through
- A viscous, expanded wet foam is pushed through the retentate tangential to the filter surface capturing the particles

The InnovaPrep Concentration Process



Frame: 0001

Wet Foam Elution

- Mode of Operation
 - “Expanded Liquid”
 - Increased Viscosity – reduced channeling (Yan et al., 2006)
 - Moves as rigid body w/ narrow ($<10\ \mu\text{m}$) lubricating layer (Briceno and Joseph, 2003; Tisne et al., 2004)
 - Bursting bubbles
- Improved Elution
 - Lower minimum elution volumes
 - Higher elution efficiencies
 - Largely unaffected by sample matrix
 - Dramatic advantage for multi-fiber & flat membrane modules
- Quickly Breaks Down into a Liquid



InnovaPrep Biological Concentrators



HCl-40 (patents pending)



HSC-40 (patents pending)

Performance

Performance Metrics

(for Autonomous Biothreat Detection Systems)

- Speed
 - Generally less than 5 to 10 minutes is desired
 - Often less than 1 or 2 minutes is desired
- Reusability
 - Reuse for hundreds of test runs is desired
 - Must withstand inhibitor materials (blinding)
- Efficiency/Concentration Factor
 - >50% efficiency is desired
 - >5x concentration factor desired
- Low instrument cost
- Low per sample cost

Speed

Speed is tied to:

- Membrane pore size
- Membrane surface area
- Sample matrix/particle loading & makeup

*For most biodefense applications bacteria, viruses, free DNA, proteins and toxins are concentrated into a single fraction. Approximates rates are:

~1 mL/min to concentrate into <200 uL

~3 mL/min to concentrate into <400 uL

Reusability

Reusability is primarily tied to the parameters listed below

- Membrane surface area
- Sample matrix/particle loading and makeup

*Reuse for up to 250 test runs has been demonstrated

Efficiency

Concentration and Recovery of Carboxylate Polystyrene Microspheres

| Particle Size, μm | Represented Particle | Average Recovery |
|------------------------------|--------------------------|------------------|
| 4.5 | Agglomerates of bacteria | 97% |
| 1.0 | Single bacteria | 90% |
| 0.05 | Viruses | 79% |
| 0.025 | DNA, Lower limit-viruses | 60% |

Efficiency

Concentration and Recovery of Microorganisms

| Species | Concentration efficiency |
|--|--------------------------|
| <i>Bacillus thuringiensis israeliensis</i> (Bti) | 97 % |
| <i>Bacillus atrophaeus</i> (Bg) | 95% |
| <i>Escherichia coli</i> (Ec) | 103% |
| MS-2 bacteriophage | 92% |
| Bg DNA | 60% |

Cost

- Instrument is a relatively simple, compact electro-mechanical instrument and thus cost is relatively low.
- Sharing of common components can reduce cost further.
- Per sample cost is very low due to reuse of concentration cells.

Integration with the Evogen DryClone

Integration for Concentration and Fluid Reuse




Integration with the Evogen DryClone


- Cyclonic Collector w/ automated wet rinse
- Concentrate Volume Range Tested - 175 to 250 μL
- Demonstrated Concentration Factor - 3.1 to 20.1
- Demonstrated Concentration Efficiency – 78% to 115%
- Fluid Reuse Demonstrated



Collection Performance vs. Extract Volume
(at 1 μm particle size)



Fluid Reuse





Integration for DHS AESaP

US Department of Homeland Security
Automated Environmental Sample Preparation Program
(AESaP)



Passage & Concentration Recoveries


- Flat membrane
- Efficiency for passage of viruses, proteins, DNA
- Efficiency for concentration of bacteria




High Yield of Viable Vegetative Bacteria

- Flat membrane
- Yield by qPCR
- Viability by plating

| | Innovaprep Yield | Viability (relative to input) | Process yield of viable cells |
|-----------|------------------|-------------------------------|-------------------------------|
| E. coli | 94 ± 15% | 77 ± 23% | 71 ± 16% |
| S. aureus | 102 ± 3% | 98 ± 3% | 100 ± 6% |



InnovaPrep input



Cell Fraction output

Integration with MFSI M-BAND



Microfluidic Systems Inc.'s (MFSI) Microfluidic Bioagent
Networked Detector (M-BAND)



Interferent Testing




Extraction Efficiency



Interferent Testing

Sample Processing Time



Contact Information

Dave Alburty

dalburty@innovaprep.com

Phone: 816.619.3375

Andy Page

apage@innovaprep.com

Phone: 816.868.6204

Zachary Packingham

zpack@innovaprep.com

Phone: 816.739.1905

InnovaPrep LLC

132 E. Main St.

Drexel, Missouri 64742

www.innovaprep.com

